

Workplace Drug & Alcohol Testing

Did you know?



**ROYAL
MEDICAL**

- Substance abuse costs Australian industry \$3.7 billion annually
- 2.5% of the Australian workforce has gone to work under the influence of illicit drugs
- 17% of the Australian workforce has used illicit drugs in the last 12 months
- Individuals who use drugs are 3.5 times more likely to injure themselves
- Approximately 7.5% of work-related injuries in Australia are alcohol related

With local statistics like these, is it any wonder that Australian industry is increasingly screening its employees for impairment and recent drug use?

URINE DRUG SCREENING



The DrugCheck® NxStep Onsite Test Cup

Complying with AS/NZS 4308:2008, this fully integrated cup test offers reliable results in just a few minutes. It comes with a resealable lid and a wider "female-friendly" opening.

ORDER CODE	CONFIGURATION	Other customised configurations are available upon request.
60601-6	AMPHETAMINES 300, BENZODIAZEPINES 200, COCAINE 300, METHAMPHETAMINES 300, OPIATES 300, THC 50 + 6 Adulterant Check + Temperature Strip.	
60702-6A	AMPHETAMINES 300, BENZODIAZEPINES 200, COCAINE 300, METHAMPHETAMINES 300, OPIATES 300, OPIATES 2000, THC 50 + Alcohol + 6 Adulterant Check + Temperature Strip.	
60930-6A	AMPHETAMINES 300, BENZODIAZEPINES 200, BUPRENORPHINE 10, COCAINE 300, METHADONE 300, METHAMPHETAMINES 300, OPIATES 300, OXYCODONE 100, THC 50 + Alcohol + 6 Adulterant Check + Temperature Strip.	



The DrugCheck® NxScan Onsite Test Cup



One of the simplest instant drug tests is also available in a flat sided cup for easy photocopying or scanning with the VeriCheck system. The VeriCheck system generates donor reports that include colour images of the actual test cup strips, along with demographics and comments entered by the testing officer.

VeriCheck®



ORDER CODE	CONFIGURATION	Other customised configurations are available upon request.
70930-6A	AMPHETAMINES 300, BENZODIAZEPINES 200, BUPRENORPHINE 10, COCAINE 300, METHADONE 300, METHAMPHETAMINES 300, OPIATES 300, OXYCODONE 100, THC 50 + Alcohol + 6 Adulterant Check + Temperature Strip.	
70601-6	AMPHETAMINES 300, BENZODIAZEPINES 200, COCAINE 300, METHAMPHETAMINES 300, OPIATES 300, THC 50 + 6 Adulterant Check + Temperature Strip.	

ORAL FLUID DRUG SCREENING



SalivaScan™

The SalivaScan™ Oral Fluid Drug Test is a sensitive screen for the detection of illicit drugs in saliva. A saturation indicator window turns red when sufficient saliva has been collected. This instant drug test has been especially optimised to detect parent THC compounds. Results are available in just a few minutes. Confirmation collection kits are included in each box of 25. Complies with AS4760.

ORDER CODE	CONFIGURATION	Other customised configurations are available upon request.
80600	AMPHETAMINES, BENZODIAZEPINES, COCAINE, METHAMPHETAMINES, OPIATES, THC + 2 CONFIRMATION TEST KITS	
80502	AMPHETAMINES, COCAINE, METHAMPHETAMINES, OPIATES, THC + 2 CONFIRMATION TEST KITS	

ALCOHOL TESTING

AlcoCheck™ Breath Alcohol

The AlcoCheck™ Breath Alcohol Test offers users a simple, convenient method to detect alcohol use and to estimate levels of alcohol intoxication using chemical alcohol detection technology with 50 years of proven success in law enforcement and corrections.

ORDER CODE: 40905



ORDER CODE: 40902

To perform a test with the AlcoCheck™ Breath Alcohol Test, the user squeezes the tube to crush the ampoule in the test, and then blows into the appropriate end of the test for 12 seconds. Yellow crystals remain unchanged if alcohol is not present at the test's threshold.

Alcohol in the breath turns yellow crystals to a pale blue or green colour within two minutes if the concentration meets or exceeds the test's threshold.

SC100 Semiconductor

- State-of-the-art semiconductor oxide sensor
- Precalibrated, easily replaceable sensor module
- Onboard memory
- Cumulative test counter
- 4-digit LCD backlight display
- Advanced blow flow detection
- Low battery indicator/Auto power-off



ORDER CODE: SC100

FC200 Fuel Cell

- Advanced European fuel cell sensor technology
- Cumulative test counter
- Onboard memory
- Portable printer capable
- 4-digit LCD backlight display
- Advanced blow flow detection
- Low battery indicator/Auto power-off



ORDER CODE: FC200

MOUTHPIECES (x50)
ORDER CODE: ATA010

K2 /SPICE ONSITE TEST NOW AVAILABLE!

Royal Medical Supplies is pleased to announce the availability of an on-site test for the detection of Synthetic Cannabinoid metabolites in urine. Initially available as a single dip-test, the test strip will quickly be adapted to the DrugCheck® NxStep and NxScan drug cups. As the test strip evolves, more of the chemical compounds commonly found in the Synthetic Cannabinoid products will be detected.

Street names

Kronic, K2, Spice, Bliss, Black Mamba, Bombay Blue, Cloud Nine, Fake Weed, Genie, Spice, Zohai.

Looks like

K2 is typically sold in small glossy cardboard sachet, containing a plastic bag of dried leaves or herbs and marketed as incense that can be smoked. It is said to resemble potpourri.

Methods of abuse

K2 products are usually smoked in joints or pipes, but some users make it into a tea. Psychological effects are similar to those of marijuana and include paranoia, panic attacks, and giddiness.



ORDER CODE: 30122-C



**ROYAL
MEDICAL**

"Improving healthcare and workplace safety through innovation".

Australian Importer and Provider of Technical Support

T: (+612) 9939 4122
F: (+612) 9939 4133
E: info@royalmedical.com.au

Our websites:

www.royalmedical.com.au
www.royalmedical.co.nz
www.drugcheck.com.au
www.drug-test.com.au